

Blood Glucose Monitoring Policy and Procedure

Charlottesville-Albemarle Rescue Squad, Inc.

Weekly Control Documentation (3 weeks per sheet)

Bayer Ascensia Contour Blood Glucose Monitoring System

Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
<i>Ex:</i> 10/9/07	<i>Ex:</i> #1	<i>Ex:</i> 3078148	<i>Ex:</i> 7AC3C53	<i>Ex:</i> 01/09	<i>Ex:</i> 100-138
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service
<i>Ex:</i> 1508	<i>Ex:</i> 03/09	<i>Ex:</i> 107	<i>Ex:</i> Pass	<i>Ex:</i> N/A	<i>Ex:</i> N/A
XXXXX	XXXXX	XXXXX	XXXXX	XXXXX	XXXXX
Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service
XXXXX	XXXXX	XXXXX	XXXXX	XXXXX	XXXXX
Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service

October 9, 2007

Revised: _____

Jim Miller, RN, CFRN, NREMT-P